

Mouse/Rat Toes (and other small animals)

<u>Station</u>	<u>Reagent</u>	<u>Time</u>	<u>Temp.</u>
1	80% alcohol	60 minutes	RT
2	95% alcohol	60 minutes	RT
3	95% alcohol	90 minutes	RT
4	95% alcohol	120 minutes	RT
5	100% alcohol	120 minutes	RT
6	100% alcohol	120 minutes	RT
7	100% alcohol	180 minutes	RT
8	Xylene	30 minutes	RT
9	Xylene	30 minutes	RT
10	Xylene	60 minutes	RT
11	Paraffin	90 minutes	60°C
12	Paraffin	120 minutes	60°C

Decalcify 1-3 days in 8% HCl/8% Formic Acid Solution