

**Mouse/Rat Skin**

| <u>Station</u> | <u>Reagent</u>                | <u>Time</u> | <u>Temp.</u> | <u>P/V</u> |
|----------------|-------------------------------|-------------|--------------|------------|
| 1              | 10% neutral buffered formalin | 60 minutes  | RT           | Yes        |
| 2              | 70% alcohol                   | 60 minutes  | RT           | Yes        |
| 3              | 80% alcohol                   | 60 minutes  | RT           | Yes        |
| 4              | 80% alcohol                   | 60 minutes  | RT           | Yes        |
| 5              | 95% alcohol                   | 60 minutes  | RT           | Yes        |
| 6              | 95% alcohol                   | 60 minutes  | RT           | Yes        |
| 7              | 100% alcohol                  | 60 minutes  | RT           | Yes        |
| 8              | 100% alcohol                  | 60 minutes  | RT           | Yes        |
| 9              | Xylene                        | 60 minutes  | RT           | Yes        |
| 10             | Xylene                        | 60 minutes  | RT           | Yes        |
| 11             | Paraffin                      | 60 minutes  | 58°C         | Yes        |
| 12             | Paraffin                      | 60 minutes  | 58°C         | Yes        |
| 13             | Paraffin                      | 60 minutes  | 58°C         | Yes        |